

SO-LOW CHEST STYLE FREEZER TO -85°C MODEL C85-5

**5 Cubic Ft. (142 liters) Chest – Temperature Range: -40°C to -85°C
Top Opening**

Specifications:

| | |
|----------------------|--|
| * Capacity | 5 cubic ft. / 142 liters |
| * Temperature Range | -40°C to -85°C |
| * Dimensions | |
| * Exterior W x D x H | 40.5" x 31" x 47" 102 x 81 x 119 (cm) |
| * Interior W x D x H | 30" x 18" x 16" 76 x 46 x 40 (cm) |

TEMPERATURE CONTROL

Microprocessor control displays set point and chamber temperature.

ALARM SYSTEM WITH RELAY

The battery operated alarm system will emit an audible and visual signal when there is a mechanical or electrical failure. The alarm has an over and under temperature setting, alarm silencing switch, and battery test switch. Also provided is a relay for remote alarm hook-up.

REFRIGERATION SYSTEM

A So-Low cascade system with two hermetic compressors. CFC & HCFC free refrigerants

CONSTRUCTION

Chamber is 14-ga. Zinc coated galvanized steel. Exterior is 16-ga. steel. Power coated cool gray.

VOLTAGE

115 volts 60 hertz, 1 phase.
15 AMP dedicated circuit required.
ETL electrical approval.

TECHNICAL SUPPORT

Supplied by Nation-Wide-Network of service companies.

Specifications subject to change without notice

