

SO-LOW CHEST STYLE FREEZER TO -85°C

MODEL C85-9

9 Cubic Ft. (255 liters) Chest – Temperature Range: -40°C to -85°C

Top Opening

Specifications:

* Capacity	9 cubic ft. / 255 liters
* Temperature Range	-40°C to -85°C
* Dimensions	
* Exterior W x D x H	57.5" x 29" x 46.5" 146 x 73 x 116 (cm)
* Interior W x D x H	45.5" x 16" x 20" 115 x 40 x 50 (cm)

TEMPERATURE CONTROL

Microprocessor control displays set point and chamber temperature.

ALARM SYSTEM WITH RELAY

The battery operated alarm system will emit an audible and visual signal when there is a mechanical or electrical failure. The alarm has an over and under temperature setting, alarm silencing switch, and battery test switch. Also provided is a relay for remote alarm hook-up.

REFRIGERATION SYSTEM

A So-Low cascade system with two hermetic compressors. CFC & HCFC free refrigerants

CONSTRUCTION

Chamber is 14-ga. Zinc coated galvanized steel. Exterior is 16-ga. steel. Power coated cool gray.

VOLTAGE

20 AMP dedicated circuit required.
ETL electrical approval.

TECHNICAL SUPPORT

Supplied by Nation-Wide-Network of service companies.

Specifications subject to change without notice

