

SO-LOW MODEL CH25-15

CHEST STYLE FREEZER TO -25°C

Capacity [Cu. Ft.]	15 Cubic Feet
Capacity [Liters]	424 Liters
Temperature Range [°C]	-0°C to -25°C
Temperature Range [°F]	+32°F to -13°F

Dimensions (W x D x H)

Exterior [Inches] 55.5" x 29.75" x 34.38"

Exterior [Centimeters] 140 x 75 x 87

Interior [Inches] 49.25" x 21.63" x 25.75"

Interior [Centimeters] 125 x 54 x 65

Electrical:

Voltage: 115 Volts 60 Hertz 1 Phase

Amperage: 15 Amp Dedicated Line

2 Running Amps

NEMA Configuration NEMA 5-15

Certification(s) ETL Listed

Construction:

Temperature Control	Digital Display
	High / Low Alarm
	Min / Max Temperature
	Dry Contacts
Alarm Relay	Manual
Defrost	Top Opening Lid
Opening	No
Lockable Lid/Door	Steel
Interior / Exterior	Polyurethane Foam
Insulation	Hermetic Compressor
Compressors	R600
Primary Refrigerant(s)	Gravity Flow
Air Circulation	4, Lockable
Casters	None
Shelving	600
Approximate BTU/Hour	55
Approximate Decibels	

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 1 Year Parts & 90 Days Labor