

## ***SO-LOW CHEST STYLE FREEZER TO -25°C***

### ***MODEL CH25-22***

**22 Cubic Ft. (708 liters) Chest – Temperature Range: 0°C to -25°C**

#### **Top Opening**

#### ***Specifications:***

* Capacity	22 cubic ft.
* Temperature Range	0°C to -25°C
* Dimensions	
* Exterior w x d x h	73.25" x 29.5" x 38.5"
* Interior w x d x h	66" x 20.5" x 27.5"

#### **CONTROL**

Digital Temperature control with two Temperature displays. One display shows the set temperature, the other shows the actual temperature in the freezer.

#### **ALARM**

Battery operated alarm system provides An audible and visual signal in a high/low temperature alarm condition. The alarm will also signal in the event of a power failure. Also provided is an alarm silencing and battery test feature.

MANUAL DEFROST

CASTERS FOR EASY MOBILITY

HERMETIC MOTOR COMPRESSOR

115 VOLT, 60 HERTZ, 1 PHASE

WARRANTY One Year Parts  
90 Days Labor – Shipping Weight 450 lbs.

#### **TECHNICAL SUPPORT**

Supplied by Nation-Wide-Network of service companies.

