

## *SO-LOW CHEST STYLE FREEZER TO -40°C*

### *MODEL CH40-5*

**5 Cubic Ft. (142 liters) Chest – Temperature Range: 0°C to -40°C**  
**Top Opening**

#### ***Specifications:***

* Capacity	5 cubic ft. / 142 liters
* Temperature Range	0°C to -40°C
* Dimensions	
* Exterior w x d x h	29" x 24" x 33.63"
* Interior w x d x h	24" x 15.25" x 26"

#### **CONTROL**

Digital Temperature control with two Temperature displays. One display shows the set temperature, the other shows the actual temperature in the freezer.

#### **ALARM**

Battery operated alarm system provides An audible and visual signal in a high/low temperature alarm condition. The alarm will also signal in the event of a power failure. Also provided is an alarm silencing and battery test feature.

MANUAL DEFROST

CASTERS FOR EASY MOBILITY

HERMETIC MOTOR COMPRESSOR

115 VOLT, 60 HERTZ, 1 PHASE

WARRANTY One Year Parts  
90 Days Labor

#### **TECHNICAL SUPPORT**

Supplied by Nation-Wide-Network of service companies.

