

SO-LOW MODEL CH45-3

MINI CHEST STYLE FREEZER TO -45°C

Capacity [Cu. Ft.]	3.5 Cubic Feet
Capacity [Liters]	99 Liters
Temperature Range [°C]	-0°C to -45°C
Temperature Range [°F]	+32°F to -49°F

Dimensions (W x D x H)

Exterior [Inches] 22" x 23.5" x 34"

Exterior [Centimeters] 55 x 59 x 86

Interior [Inches] 15" x 17" x 25.87"

Interior [Centimeters] 40 x 43 x 65

Electrical:

Voltage: 115 Volts 60 Hertz 1 Phase

Amperage 15 Amp Dedicated Line

2 Running Amps

NEMA Configuration NEMA 5-15

Certification(s) ETL Listed

Construction:

Temperature Control	Digital Display
	High / Low Alarm
	Min / Max Temperature
	Dry Contacts
Alarm Relay	Manual
Defrost	Top Opening Lid
Opening	No
Lockable Lid/Door	Steel
Interior / Exterior	Polyurethane Foam
Insulation	Hermetic Compressor
Compressors	R290
Primary Refrigerant(s)	Gravity Flow
Air Circulation	4, Lockable
Casters	None
Shelving	600
Approximate BTU/Hour	55
Approximate Decibels	

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 1 Year Parts & 90 Days Labor