

SO-LOW LAB & PHARMACY REFRIGERATOR/FREEZER MODEL DH-10RFDA

Refrigerator 7 cu. ft. - Freezer 3 cu. ft.

Upright – Temperature Range: 2°C/8°C (R) / -10°C/-25°C (F) Two Doors

Specifications:

* Capacity	10 cubic ft.
* Temperature Range	2°C to 8°C / -10°C to -25°C
* Dimensions	
* Exterior w x d x h	23.75" x 25.5" x 66.75" 60 x 65 x 170 (cm)
* Fridge Interior w x d x h	20.25" x 18.5" x 28.75"
* Freezer Interior w x d x h	18" x 17.5" x 14"

Features:

Temperature Control	Two - Digital display
Temperature Alarm	Hi/Low temperature alarm With audible, and visual alarm Alarm relay dry contacts
Doors	Freezer solid, locking Refrigerator glass, locking
Casters	Four
Shelves	Three adjustable / Refrigerator Two adjustable / Freezer
Probe Port	3/8" access port hole, EA chamber
Insulation	Urethane foam
Refrigerant:	CFC free
Compressor	Two
Refrigerator Defrost	Cycle, automatic
Freezer Defrost	Manual
Voltage: Domestic	115/60/1
Supply Plug	5-15P NEMA
Shipping Weight	395 lbs.
Warranty-(Domestic)	Two year parts and labor Five year compressor warranty



Specifications subject to change without notice