

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

## **SO-LOW LAB AND PHARMACY REFRIGERATOR FREEZER COMBINATION MODEL DH-20RF**

**Refrigerator 9.1 cu. ft. -Freezer 9.1 cu. ft.**

**Upright – Temperature Range: 8°C (R) / -25°C (F) Two Doors**

### **Specifications:**

* Total Capacity		18.2 cubic ft.
* Temperature (Fridge)		2°C to 8°C
* Temperature (Freezer)		0°C to -25°C
<b>*When ordering please specify Celsius or Fahrenheit Setting*</b>		
* Dimensions		
* Exterior w x d x h		27.5" x 34" x 81.5" 70 x 86 x 207 (cm)
* Fridge Interior w x d x h		28 x 28 x 20"
* Freezer Interior w x d x h		28 x 28 x 20"

### **Features:**

Control	Digital Display with Hi/Low Alarm Alarm Relay Dry Contacts For each compartment
Defrost	Auto
Doors	Two Solid, Locking
Casters	Four, Lockable
Interior Finish	Stainless Steel
Exterior Finish	Stainless Steel
Shelves	Epoxy Coated / 2 Per Compartment
Insulation:	Polyurethane
Evaporator	Fin and tube
Compressors	1/4 hp (fridge), 1/3 hp (freezer)
Voltage:	115/60/1
Running Amps	4.5 (fridge), 7.5 (freezer)
Supply Plugs	(2) NEMA 5-15P UL listed
Shipping Weight	550 lbs.
Warranty	Three year parts and labor Compressors additional 2 years, no labor.



Specifications subject to change without notice