

SO-LOW LAB AND PHARMACY REFRIGERATOR FREEZER COMBINATION MODEL DH-48RF

Refrigerator 24 cu. ft. -Freezer 24 cu. ft.

Upright – Temperature Range: 8°C (R) / -25°C (F) Two Doors

Specifications:

- * Total Capacity 48 cubic ft.
- * Temperature (Fridge) 2°C to 8°C
- * Temperature (Freezer) 0°C to -25°C
- *When ordering please specify Celsius or Fahrenheit Setting***

*** Dimensions**

- * Exterior w x d x h 52 x 34.75 x 81.5”
132 x 88 x 207”

- * Fridge Interior w x d x h 22 x 28 x 60”
- * Freezer Interior w x d x h 22 x 28 x 60”

Features:

- | | |
|----------------------|--|
| Control | Digital Display with Hi/Low Alarm
Alarm Relay Dry Contacts
Min/Max Memory
Power Failure Relay
For each compartment |
| Defrost | Auto |
| Doors | Two Solid, Locking |
| Casters | Four, Lockable |
| Interior Finish | Stainless Steel |
| Exterior Finish | Stainless Steel |
| Shelves | Epoxy Coated / Three Per Compartment |
| Insulation: | Polyurethane |
| Evaporator | Fin and tube |
| Compressor | 1/4 hp - 1/2 hp |
| Voltage: | 115/60/1 |
| Running Amps/Breaker | 5.7/15 – 6.0/15 |
| Supply Plug | NEMA 5-15P
UL listed |
| Shipping Weight | 625 lbs. |
| Warranty-(Domestic) | Three year parts and labor
Compressor additional
2 years, no labor. |

