

## SO-LOW LAB AND PHARMACY REFRIGERATOR FREEZER COMBINATION MODEL DH-70RF

**Refrigerator 46.6 cu. ft. -Freezer 24 cu. ft.**

**Upright – Temperature Range: 4°C (R) / -20°C (F) Three Doors**

### Specifications:

- \* Total Capacity 70 cubic ft.
- \* Temperature (Fridge) 4°C
- \* Temperature (Freezer) -20°C
- \*When ordering please specify Celsius or Fahrenheit Setting\***

### \* Dimensions

- \* Exterior w x d x h 78 x 34.75 x 81.5”  
198 x 88 x 207 (cm)

- \* Fridge Interior w x d x h 48 x 28 x 60”
- \* Freezer Interior w x d x h 28 x 28 x 20”

### Features:

- |                      |  |
|----------------------|--|
| Control              | Digital Display with Hi/Low Alarm<br>Alarm Relay Dry Contacts<br>Min/Max Memory<br>Power Failure Relay<br>For each compartment |
| Defrost              | Auto   |
| Doors                | Three Solid, Locking   |
| Casters              | Four, Lockable   |
| Interior Finish      | Stainless Steel  |
| Exterior Finish      | Stainless Steel  |
| Shelves              | Epoxy Coated / 3 Per Compartment   |
| Insulation:          | Polyurethane   |
| Evaporator           | Fin and tube   |
| Compressor           | 1/3 hp - 1/2 hp  |
| Voltage:             | 115/60/1   |
| Running Amps/Breaker | 9.4/15 – 10.6/15   |
| Supply Plug          | NEMA 5-15P<br>UL listed  |
| Shipping Weight      | 950 lbs.   |
| Warranty-(Domestic)  | Three year parts and labor<br>Compressor additional<br>2 years, no labor.  |

