

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

SO-LOW LAB & PHARMACY FREEZER *MODEL DH29-74SD*

“SETTING TEMPERATURE STANDARDS IN EVERY LAB”

74 Cubic Ft. (2095 liters) Upright

Temperature Range: -18°C to -25°C / 0°F to -13°F
3 Solid Doors

Specifications:

* Capacity 74 cubic ft. / 2095 liters
* Temperature Range -18°C to -25°C / 0°F to -13°F
When ordering please specify Celsius or Fahrenheit Setting

* Dimensions
* Exterior w x d x h 78" x 34" x 81.5"
198 x 86 x 207 (cm)

* Interior w x d x h 74" x 28" x 60"
188 x 71 x 152 (cm)

Features:

| | |
|---------------------------|---|
| Temperature Control Keyed | Digital display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay |
| Access Port | 1/2" Diameter Top Back Wall |
| Door | Three doors solid – lockable |
| Casters | Heavy duty, swivel |
| Shelves | 9 Epoxy coated adjustable |
| Interior Lights | Two 40W |
| Insulation | Polyurethane CFC Free |
| Refrigerant: | R-404A |
| Compressor | 1 hp |
| Air Circulation | Forced air |
| Evaporator | Fin and tube, forced air |
| Thermostat | Adjustable |
| Defrost | Timed automatic |
| Voltage: | 230/115/60/1 |
| Running Amps/Breaker | 8.4/20 |
| Supply Plug | L14-20 NEMA |
| Shipping Weight | 900 lbs. |
| Warranty-(Domestic) | Three year parts and labor Compressor additional 2 years, no labor. |



Specifications subject to change without notice.