

SO-LOW LAB & PHARMACY FREEZER MODEL DH30-49SDSS

49 Cubic Ft. Upright – Temperature Range: -18°C to -25°C / 0°F to -13°F

Two Solid Doors

Specifications:

- * Capacity 49 cubic ft. / 1387 liters
- * Temperature Range -18°C to -25°C / 0°F to -13° F
- *When ordering please specify Celsius or Fahrenheit Setting***
- * Dimensions
- * Exterior w x d x h 52" x 33.75" x 81"
132 x 85 x 206 (cm)
- * Interior w x d x h 47.5" x 28.5" x 55"
121 x 72 x 140 (cm)

Features:

- | | |
|----------------------|---|
| Temperature Control | Digital display
W/ Hi-Low Alarm
Alarm Relay Dry Contacts
Min/Max Memory
Power Failure Relay |
| Access Port | 1/2" Diameter Top Back Wall |
| Doors | Two solid locking |
| Interior Finish | Stainless Steel |
| Exterior Finish | Stainless Steel |
| Casters | Four, Locking |
| Shelves | Six |
| Interior Lighting | One 40-watt fluorescent |
| Insulation | Urethane foam |
| Refrigerant: | R404A |
| Compressor | 3/4 hp |
| Air Circulation | Forced air |
| Evaporator | Fin and tube, forced air |
| Defrost | Automatic |
| Voltage: Domestic | 115/60/1 |
| Running Amps/Breaker | 13.8/20 |
| Supply Plug | 5-20P NEMA |
| Shipping Weight | 610 lbs. |
| Warranty-(Domestic) | Three year parts and labor
Compressor additional
2 years, no labor. |

