

MODEL DH4-12GD

LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.]	9.19 Cu. Ft.
Capacity [Liters]	260 Liters
Temperature Range [°C]	2°C to 8°C
Temperature Range [°F]	36°F to 46°F

Dimensions	Width x Depth x Height
Exterior [Inches]	25" x 22.75" x 63.25"
Interior [Inches]	22" x 16.5" x 43.75"
Exterior [Centimeters]	63 x 57 x 160
Interior [Centimeters]	55 x 41 x 111

Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Cycle
Refrigerant	R290
Exterior	Painted Steel
Interior	Painted Steel
Insulation	Polyurethane Foam
Lockable Door	1 Door, Keyed Lock
Door Material	Glass
Interior Light	Adjustable LED
Porthole	1/2" Diameter (Top Back Wall)
Casters	4
Storage	4 Adjustable Shelves
Shipping Weight	400 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 4.0A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed Energy Star Rated

Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Display and Reset
Alarm Relay	Dry Contacts

Warranty 3 Years Parts & Labor.
Additional 4 Years on Compressor Parts, No Labor.

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

