

## MODEL DH4-15GD

### LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

<b>Capacity [ Cu. Ft. ]</b>	<b>15.4 Cu. Ft.</b>
<b>Capacity [ Liters ]</b>	<b>436 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>2°C to 8°C</b>
<b>Temperature Range [ °F ]</b>	<b>36°F to 46°F</b>
Dimensions	Width x Depth x Height
Exterior [ Inches ]	29.5" x 26" x 78.5"
Interior [ Inches ]	25.5" x 19.25" x 54.25"
Exterior [ Centimeters ]	74 x 66 x 199
Interior [ Centimeters ]	64 x 48 x 137
Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Cycle
Refrigerant	R290
Exterior	Painted Steel
Interior	Painted Steel
Insulation	Polyurethane Foam
Lockable Door	1 Door, Keyed Lock
Door Material	Glass
Interior Light	Adjustable LED
Porthole	1/2" Diameter ( Top Back Wall )
Casters	4
Storage	5 Adjustable Shelves
Shipping Weight	440 lbs.

#### Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line   Running	15 Amp Dedicated   7.0A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed Energy Star Rated

#### Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Display and Reset
Alarm Relay	Dry Contacts

#### Warranty

3 Years Parts & Labor.

Additional 4 Years on Compressor Parts, No Labor.

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

