

SO-LOW LABORATORY REFRIGERATOR MODEL DH4-15SGD

“SETTING TEMPERATURE STANDARDS IN EVERY LAB”

15 Cubic Ft. Upright – Temperature Range: 2°C to 8°C

Two Glass Sliding Doors

Specifications:

* Capacity	15 cubic ft.
* Temperature Range	2°C to 8°C
* Dimensions	
* Exterior w x d x h	36" x 24" x 65" 91 x 60 x 165 (cm)
* Interior w x d x h	33 x 21" x 41" 83 x 53 x 104 (cm)

Features:

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts
Access Port	1/2" Diameter Top Back Wall
Exterior	White coated steel
Interior	White coated steel
Doors	Double glass – sliding Positive seal, self closing
Casters	Four, locking
Shelves	Three
Interior Lighting	One fluorescent light
Insulation	Foamed in place polyurethane
Refrigerant:	R134A
Compressor	1/3 hp
Evaporator	Fin and tube, forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps	5.8/15
Supply Plug	5-15P NEMA
Shipping Weight	450 lbs.
Warranty-(Domestic)	Three year parts and labor Compressors additional 2 years, no labor.

