

## ***SO-LOW CHROMATOGRAPHY REFRIGERATOR MODEL DH4-23GDCHSS***

**23 Cubic Ft. Upright – Temperature Range: 2°C to 8°C**

**Glass Door**

### ***Specifications:***

* Capacity	23 cubic ft.
* Temperature Range	2°C to 8°C
* Dimensions	
* Exterior w x d x h	27" x 33.75" x 81" 68 x 85 x 206 (cm)
* Interior w x d x h	24" x 28.5" x 55" 61 x 72 x 140 (cm)

### ***Features:***

Temperature Control	Digital display W/ Hi-Low Alarm Alarm Relay Dry Contacts
Interior Outlet Access Ports	115V Duplex 5A LH & RH Ends of cabinet 2" diameter
Doors	Single glass, locking
Interior Finish	Stainless Steel
Exterior Finish	Stainless Steel
Casters	Four, locking
Shelves	Three
Interior Lighting	One fluorescent
Insulation	Polyurethane foam
Refrigerant:	R134A CFC free
Compressor	1/3 hp
Condenser	Fin and tube
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	5.8/15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	440 lbs.
Warranty-(Domestic)	Three year parts and labor Compressors additional 2 years, no labor.

