

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

## *SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DH4-23SD*

**23 Cubic Ft. Upright – Temperature Range: 2°C to 8°C**

**Solid Door**

### *Specifications:*

* Capacity	23 cubic ft.
* Temperature Range	2°C to 8°C
* Dimensions	
* Exterior w x d x h	27" x 34" x 81" 68 x 86 x 206 (cm)
* Interior w x d x h	24" x 28.5" x 55" 61 x 72 x 140 (cm)

### *Features:*

Temperature Control	Digital display W/ Hi-Low Alarm Alarm Relay Dry Contacts
Access Port	1/2" Diameter Top Back Wall
Doors	Single Solid, locking
Interior finish	White powder coat
Casters	Four, lockable
Shelves	Five
Interior Lighting	One 20-watt fluorescent
Insulation	Polyurethane foam
Refrigerant:	R134A CFC free
Compressor	1/3 hp
Condenser	Fin and tube
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Thermostat	Adjustable
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	5.8/15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	440 lbs.
Warranty-(Domestic)	Three year parts and labor Compressors additional 2 years, no labor.



Specifications subject to change without notice