

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

## *SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DH4-27GD*

**27 Cubic Ft. (764 liters) Upright – Temperature Range: 2°C to 8°C  
Glass Door**

### *Specifications:*

* Capacity	27 cubic ft. / 764 liters
* Temperature Range	2°C to 8°C
* Dimensions	
* Exterior w x d x h	30" x 32" x 81" 76 x 81 x 206 (cm)
* Interior w x d x h	27" x 28.5" x 55" 69 x 72 x 140 (cm)

### *Features:*

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts
Access Port	1/2" Diameter Top Back Wall
Door	Single glass, locking
Interior finish	White coated steel
Casters	Four, locking
Shelves	Five
Interior Lighting	One 20-watt fluorescent
Insulation	Polyurethane foam
Refrigerant:	R134A CFC free
Compressor	1/3 hp
Condenser	Fin and tube
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	5.8/15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	440 lbs.
Warranty-(Domestic)	Three year parts and labor Compressors additional 2 years, no labor.



Specifications subject to change without notice