

## MODEL DH4-30SGD

### LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

<b>Capacity [ Cu. Ft. ]</b>	<b>30 Cu. Ft.</b>
<b>Capacity [ Liters ]</b>	<b>849 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>2°C to 8°C</b>
<b>Temperature Range [ °F ]</b>	<b>36°F to 46°F</b>

Dimensions	Width x Depth x Height
Exterior [ Inches ]	47.25" x 30.25" x 81"
Interior [ Inches ]	43" x 22" x 54"
Exterior [ Centimeters ]	120 x 76 x 205
Interior [ Centimeters ]	109 x 55 x 137

Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Cycle
Refrigerant	R290
Exterior	Painted Steel
Interior	Painted Steel
Insulation	Polyurethane Foam
Door	2 Sliding Doors

<b>Accessibility</b>	<b>Doors are Not Lockable</b>
Door Material	Glass
Interior Light	Adjustable LED
Porthole	1/2" Diameter
Casters	4, Locking
Storage	10 Adjustable Shelves
Shipping Weight	500 lbs.

#### Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line   Running	15 Amp Dedicated   8.3A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed Energy Star Rated

#### Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Display and Reset
Alarm Relay	Dry Contacts

**Warranty** 3 Years Parts & Labor.  
Additional 4 Years on Compressor Parts, No Labor.

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

