

SO-LOW CHROMOTOGRAPHY REFRIGERATOR MODEL DH4-45SGDCH

45 Cubic Ft. (1286 liters) Upright – Temperature Range: 2°C to 8°C

Two Glass Sliding Doors

Specifications:

* Capacity	45 cubic ft. / 1286 liters
* Temperature Range	2°C to 8°C
* Dimensions	
* Exterior w x d x h	52" x 32" x 81" 132 x 81 x 206 (cm)
* Interior w x d x h	49" x 27" x 55" 125 x 69 x 140 (cm)

Features:

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts
Interior	White powder-coated steel
Interior Outlet Ports	115V Duplex 5A Access ports on LH & RH Ends of cabinet – 2" diameter
Doors	Double glass – sliding Positive seal, self closing Triple pane thermal glass
Shelves	Ten, epoxy coated wire
Interior Lighting	One fluorescent light
Insulation	Polyurethane
Refrigerant:	R134A
Compressor	1/3 hp
Evaporator	Fin and tube, forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	8.8/15
Supply Plug	5-15P NEMA
Shipping Weight	610 lbs.
Warranty-(Domestic)	Two year parts and labor Compressors additional 3 years, no labor.

