

## MODEL DH4-49GD

### LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

**Capacity [ Cu. Ft. ]** 49 Cu. Ft.  
**Capacity [ Liters ]** 1387 Liters  
**Temperature Range [ °C ]** 2°C to 8°C  
**Temperature Range [ °F ]** 36°F to 46°F

Dimensions  
Exterior [ Inches ] 52" x 32" x 81"  
Interior [ Inches ] 47.5" x 28.5" x 55"  
Exterior [ Centimeters ] 132 x 81 x 206  
Interior [ Centimeters ] 120 x 72 x 140

Compressor Hermetic Compressor  
Air Circulation Forced Air  
Defrost Method Cycle  
Refrigerant R290  
Exterior Painted Steel  
Interior Painted Steel  
Insulation Polyurethane Foam  
Lockable Doors Two Doors, Keyed Lock  
Door Material Glass  
Interior Light Adjustable LED  
Porthole 1/2" Diameter ( Top Back Wall )  
Casters 4, Locking  
Storage 10 Adjustable Shelves  
Shipping Weight 610 lbs.

#### Electrical

Voltage: 115 Volts 60 Hertz 1 Phase  
Amperage Line | Running 15 Amp Dedicated | 5.0A  
NEMA Configuration NEMA 5-15  
Certification(s) UL Listed  
Energy Star Rated

#### Controller

Display Type Microprocessor, Digital  
Alarms High/Low Temp  
Min / Max Temperature Display and Reset  
Alarm Relay Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



**Warranty** 7 Years Parts, Labor & Compressor Part Warranty.