

SO-LOW CHROMATOGRAPHY REFRIGERATOR MODEL DH4-49GDCH-SS

49 Cubic Ft. (1387 liters) Upright – Temperature Range: 2°C to 8°C

Two Swinging Glass Doors

Specifications:

* Capacity	49 cubic ft. / 1387 liters
* Temperature Range	2°C to 8°C
* Dimensions	
* Exterior w x d x h	52" x 33.75" x 81" 132 x 85 x 206 (cm)
* Interior w x d x h	47.5" x 28.5" x 55" 120 x 72 x 140 (cm)

Features:

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts
Interior Outlet Ports	115V Duplex 5A Access ports on LH & RH Ends of cabinet – 2" diameter
Doors	Two hinged glass, locking
Interior Finish	Stainless Steel
Exterior Finish	Stainless Steel
Casters	Four, Locking
Shelves	Six
Interior Lighting	One 40-watt fluorescent
Insulation	Polyurethane foam
Refrigerant:	R134A
Compressor	1/3 hp
Condenser	Fin and tube
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	8.8/15
Supply Plug	5-15P NEMA
Shipping Weight	610 lbs.
Warranty-(Domestic)	Two year parts and labor Compressors additional 3 years, no labor.



Specifications subject to change without notice