

## MODEL DH4-50PT

### PASS-THRU LAB REFRIGERATOR 2°C TO 8°C

<b>Capacity [ Cu. Ft. ]</b>	<b>50 Cu. Ft.</b>
<b>Capacity [ Liters ]</b>	<b>1415 Liters</b>
<b>Temperature Range [ °C ]</b>	<b>2°C to 8°C</b>
<b>Temperature Range [ °F ]</b>	<b>36°F to 46°F</b>
Dimensions	Width x Depth x Height
Exterior [ Inches ]	52" x 37.5" x 81.5"
Interior [ Inches ]	48" x 30" x 59"
Exterior [ Centimeters ]	132 x 95 x 207
Interior [ Centimeters ]	121 x 76 x 149
Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Cycle
Refrigerant	R290
Exterior	Stainless Steel
Interior	Stainless Steel
Insulation	Polyurethane Foam
Lockable Door	4 Doors, Keyed Lock
Door Material	Glass ( Front & Back Doors ) Option for Solid Door(s)
Interior Light	Adjustable LED
Porthole	1/2" Diameter
Casters	4, Locking
Storage	6 Adjustable Shelves
Shipping Weight	650 lbs.

#### Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line   Running	15 Amp Dedicated   8.0A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed

#### Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Display and Reset
Alarm Relay	Dry Contacts

**Warranty** 3 Years Parts & Labor.  
Additional 4 Years on Compressor Parts, No Labor.

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

