

MODEL DH4-66SGD

LAB AND PHARMACY REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.] 66 Cu. Ft.
Capacity [Liters] 1868 Liters
Temperature Range [°C] 2°C to 8°C
Temperature Range [°F] 36°F to 46°F

Dimensions
Exterior [Inches] 75" x 32" x 81"
Interior [Inches] 72" x 27" x 55"
Exterior [Centimeters] 190 x 81 x 206
Interior [Centimeters] 182 x 69 x 140

Compressor Hermetic Compressor
Air Circulation Forced Air
Defrost Method Cycle
Refrigerant R290
Exterior Painted Steel
Interior Painted Steel
Insulation Polyurethane Foam
Lockable Door 3 Sliding Doors, Keyed Lock
Door Material Glass
Interior Light Adjustable LED
Porthole 1/2" Diameter
Casters 4, Locking
Storage 15 Adjustable Shelves
Shipping Weight 850 lbs.

Electrical

Voltage: 115 Volts 60 Hertz 1 Phase
Amperage Line | Running 15 Amp Dedicated | 6.2A
NEMA Configuration NEMA 5-15
Certification(s) UL Listed

Controller

Display Type Microprocessor, Digital
Alarms High/Low Temp
Min / Max Temperature Display and Reset
Alarm Relay Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 7 Years Parts, Labor & Compressor Part Warranty.