

SO-LOW LABORATORY REFRIGERATOR MODEL DH4-72SDSS

**72 Cubic Ft. (2095 liters) Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F
Three Hinged Solid Doors**

Specifications:

* Capacity	72 cubic ft. / 2095 liters
* Temperature Range	2°C to 8°C / 36°F to 46°F
When ordering please specify Celsius or Fahrenheit Setting	
* Dimensions	
* Exterior w x d x h	75" x 33.75" x 81" 190 x 85 x 206 (cm)
* Interior w x d x h	72" x 27" x 61" 183 x 69 x 155 (cm)

Features:

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay
Access Port	1/2" Diameter Top Back Wall
Exterior	Stainless Steel
Interior	Stainless Steel
Doors	Three solid – hinged – lockable
Shelves	Nine, Adjustable
Casters	Four, Locking
Interior Lighting	LED light
Insulation	Polyurethane CFC Free
Refrigerant	R290A CFC free
Compressor	1/2 hp
Air Circulation	Forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	7.2 /15
Supply Plug	NEMA 5-15P
Shipping Weight	850 lbs.
Warranty-(Domestic)	Three year parts and labor Compressor additional 2 years, no labor.

