

MODEL DHH-30RFX

EXPLOSION PROOF REFRIGERATOR / FREEZER COMBINATION UNIT

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

Total Capacity [Cu. Ft.] 30 Cu. Ft. (Total)
Capacity [Liters] 849 Liters (Total)

	REFRIGERATOR	FREEZER
Temperature Range [°C]	1°C to 10°C	-18°C to -25°C
Temperature Range [°F]	33°F to 55°F	0°F to -13°F

Dimensions	Width x Depth x Height	
Interior [Inches]	32" x 22.5" x 50.25"	29" x 36.25" x 77.5"
Interior [Centimeters]	81 x 57 x 127	73 x 56 x 30
Exterior [Inches]	35.75" x 36.25" x 77.5"	
Exterior [Centimeters]	90 x 92 x 196	

Compressor	1 Each, Hermetic Compressor
Refrigerator Air Circulation	Gravity Flow
Refrigerator Defrost Method	Manual
Refrigerator Storage	3 Adjustable Shelves
Freezer Air Circulation	Gravity Flow
Freezer Defrost Method	Manual
Freezer Storage	N/A
Refrigerant	R404A
Exterior	Painted Steel
Interior	Painted Steel
Insulation	Polyurethane Foam
Lockable Door	1 Door, Keyed Lock
Leveling Legs	4
Shipping Weight	400 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 6A
Electrical Plug Certification(s)	Not Supplied (Hardwired) UL Listed

Controller

Display Type	Mechanical Thermostat
--------------	-----------------------



» Meets NFPA / OSHA Guidelines for Hazardous Locations
 Class 1, Division II, Group C & D

Warranty 1 Year Parts & Labor