

## MODEL DHH-9RFFMS

# FLAMMABLE MATERIAL STORAGE REFRIGERATOR / FREEZER COMBINATION UNIT

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES

**Total Capacity [ Cu. Ft. ]** 9 Cu. Ft. (Total)  
**Capacity [ Liters ]** 255 Liters (Total)

	REFRIGERATOR	FREEZER
<b>Temperature Range [ °C ]</b>	1°C to 10°C	-15°C to -25°C
<b>Temperature Range [ °F ]</b>	33°F to 55°F	5°F to -13°F

Dimensions	Width x Depth x Height	
Interior [ Inches ]	20.25" x 17" x 27"	19" x 16" x 26"
Interior [ Centimeters ]	51 x 43 x 68	48 x 40 x 66
Exterior [ Inches ]	23.75" x 24.25" x 65.5"	
Exterior [ Centimeters ]	60 x 61 x 168	

Compressor	1 Each, Hermetic Compressor
Refrigerator Air Circulation	Forced Air
Refrigerator Defrost Method	Automatic
Refrigerator Storage	3 Adjustable Shelves

Freezer Air Circulation	Gravity Flow
Freezer Defrost Method	Manual
Freezer Storage	2 Fixed Shelves

Refrigerant	R600a
Exterior	Painted Steel
Interior	Painted Steel
Insulation	Urethane Foam
Lockable Door	1 Door Each, Keyed Lock
Leveling Legs	4
Shipping Weight	205 lbs.

### Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line   Running	15 Amp Dedicated   2A
NEMA Configuration	NEMA 5-15

### Controller

Display Type	1 Each, Microprocessor, Digital
Alarms	High/Low Temp
Min / Max Temperature	Display and Reset
Alarm Relay	Dry Contacts



» Meets NFPA / OSHA Guidelines 45 & 70 for Flammable Material Storage in General Laboratory Work Areas.

**Warranty** 2 Year Parts & Labor  
Additional 5 Years on Compressor Parts, No Labor.