

## ***SO-LOW FLAMMABLE MATERIAL STORAGE FREEZER MODEL DHH20-17SDFMS MANUAL DEFROST***

*Meets NFPA and OSHA guidelines for flammable material storage in general laboratory work areas*

### **17 Cubic Ft. Upright -25°C / Solid Door**

#### **Specifications:**

* Capacity	17 cubic ft. / 481 liters
* Temperature	-18°C to -25°C
* Dimensions	
* Exterior w x d x h	34" x 33" x 67 5/8 " 86.36 x 83.82 x 171.78 (cm)
* Interior w x d x h	27" x 17.5" x 51.5" 68.58 x 44.45 x 130.81 (cm)

#### **Features:**

*Thermostat	Exterior Adjustable Mechanical
*Doors	Single solid – key lock
*Legs	Four leveling
*Shelves	3 fixed
*Interior Lighting	None
*Insulation:	Polyurethane foam
*Refrigerant:	R-134A
*Compressor	1/3 hp
*Air Circulation	Gravity
*Evaporator	Hot wall
*Defrost	Manual
*Voltage	115-60-1
*Running Amps	5
*Supply Plug	5-15 NEMA
*Shipping weight	261 lbs.
*Approvals	UL Listed
*Warranty-(Domestic)	One year parts and labor.

