

SO-LOW FLAMMABLE MATERIAL STORAGE FREEZER MODEL DHH20-20SDFMS MANUAL DEFROST

Meets NFPA and OSHA guidelines for flammable material storage in general laboratory work areas

20 Cubic Ft. Upright -25°C / Solid Door

Specifications:

* Capacity	20 cubic ft. / 566 liters
* Temperature	-18°C to -25°C
* Dimensions	
* Exterior w x d x h	34" x 35 ¼" x 73 ¼" 86.36 x 88.9 x 185.42 (cm)
* Interior w x d x h	27" x 17.5" x 51.5" 68.58 x 44.45 x 130.81 (cm)

Features:

*Thermostat	Exterior Adjustable Mechanical
*Doors	Single solid – key lock
*Legs	Four leveling
*Shelves	4 fixed
*Interior Lighting	None
*Insulation:	Polyurethane foam
*Refrigerant:	R-134A
*Compressor	1/4 hp
*Air Circulation	Gravity
*Evaporator	Hot wall
*Defrost	Manual
*Voltage	120 VAC
*Running Amps	5
*Supply Plug	5-15 NEMA
*Shipping weight	313 lbs.
*Approvals	UL Listed
*Warranty-(Domestic)	One year parts and labor.

