

SO-LOW FLAMMABLE MATERIAL STORAGE REFRIGERATOR MODEL DHH4-17SDFMS MANUAL DEFROST

Meets NFPA and OSHA guidelines for flammable material storage in general laboratory work areas

17 Cubic Ft. Upright +4°C / Solid Door

Specifications:

* Capacity	17 cubic ft. / 481 liters
* Temperature	1°C to 10°C
* Dimensions	
* Exterior w x d x h	34" x 33" x 67 5/8" 86.36 x 83.82 x 171.75 (cm)
* Interior w x d x h	26.75" x 17.75" x 51.25" 67.95 x 45.09 x 130.18 (cm)

Features:

*Thermostat	Exterior Adjustable Mechanical
*Doors	Single solid – key lock
*Legs	Four leveling
*Shelves	1 adjustable / 2 fixed
*Interior Lighting	None
*Insulation:	Polyurethane foam
*Refrigerant:	R-134A
*Compressor	1/4 hp
*Condenser	Hot wall
*Air Circulation	Gravity
*Evaporator	Cold wall
*Defrost	Manual
*Voltage	115-60-1
*Running Amps	5
*Supply Plug	5-15 NEMA
*Shipping weight	264 lbs.
*Approvals	UL Listed
*Warranty-(Domestic)	One year parts and labor.

