

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

## *SO-LOW FLAMMABLE MATERIAL STORAGE REFRIGERATOR MODEL DHH4-20SDFMS MANUAL DEFROST*

*Meets NFPA and OSHA guidelines for flammable material storage in general laboratory work areas*

### **20 Cubic Ft. Upright +4°C / Solid Door**

#### **Specifications:**

* Capacity	20 cubic ft. / 566 liters
* Temperature	1°C to 10°C
* Dimensions	
* Exterior w x d x h	34" x 35 1/4" x 73 1/4" 86.36 x 89.54 x 186.06 (cm)
* Interior w x d x h	26.75" x 19.75" x 56.75" 67.95 x 50.17 x 144.15 (cm)

#### **Features:**

*Thermostat	Exterior Adjustable Mechanical
*Doors	Single solid – key lock
*Legs	Four leveling
*Shelves	1 adjustable / 3 fixed
*Interior Lighting	None
*Insulation:	Polyurethane foam
*Refrigerant:	R-134A
*Compressor	1/3 hp
*Condenser	Hot wall
*Air Circulation	Gravity
*Evaporator	Cold wall
*Defrost	Manual
*Voltage	120 VAC
*Running Amps	5
*Supply Plug	5-15 NEMA
*Shipping weight	316 lbs.
*Approvals	UL Listed
*Warranty-(Domestic)	One year parts and labor.



Specifications subject to change without notice