

SO-LOW FLAMMABLE MATERIAL STORAGE REFRIGERATOR MODEL DHH4-30SDFMS MANUAL DEFROST

Meets NFPA and OSHA guidelines for flammable material storage in general laboratory work areas

30 Cubic Ft. Upright +4°C / Solid Door

Specifications:

* Capacity	30 cubic ft. / 849 liters
* Temperature	1°C to 10°C
* Dimensions	
* Exterior w x d x h	35 ¾" x 35 ½" x 72 ½" 90.81 x 90.17 x 184.15 (cm)
* Interior w x d x h	32" x 22.5" x 65" 81.28 x 57.15 x 165.1 (cm)

Features:

*Thermostat	Exterior Adjustable Mechanical
*Doors	Single solid – key lock
*Legs	Four leveling
*Shelves	5 adjustable
*Interior Lighting	None
*Insulation:	Polyurethane foam
*Refrigerant:	R-134A
*Compressor	1/3 hp
*Evaporator	Hot wall
*Air Circulation	Gravity
*Defrost	Manual
*Voltage	115/60/1
*Running Amps	6
*Supply Plug	5-15 NEMA
*Shipping weight	381 lbs.
*Approvals	UL Listed
*Warranty-(Domestic)	One year parts and labor.

