

## SO-LOW LAB & PHARMACY FREEZER

### MODEL DHK30-49SD

49 Cubic Ft. Upright – Temperature Range: -15°C to -30°C

#### Solid Door

#### Specifications:

* Capacity	49 cubic ft.
* Temperature Range	-15°C to -30°C
<b>*When ordering please specify Celsius or Fahrenheit Setting*</b>	
* Dimensions	
* Exterior w x d x h	55" x 32" x 83" 137 x 81 x 212 (cm)
* Interior w x d x h	49.4" x 26.8" x 50.8" 125 x 68 x 129 (cm)

#### Features:

Temperature Control	Digital display With Hi/Low Alarm
Door	Double door, solid, locking
Access Port	1/2" Diameter
Interior Finish	Stainless
Exterior Finish	Stainless
Casters	Four Casters
Shelves	Six Adjustable
Interior Lighting	LED Light
Insulation	Polyurethane foam
Refrigerant:	R600 CFC free
Air Circulation	Forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	9 / 15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	600 lbs.
Warranty	2 Years Parts and Labor. Additional 3 years on compressor parts.

