

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DHK4-30SD-SS

30 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Solid Door

Specifications:

* Capacity	30 cubic ft.
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	33.85" x 37.80" x 81.12"
* Interior w x d x h	27.50" x 30.80" x 61.84"

Features:

Temperature Control	Digital display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay Battery back up
Door	Single solid, locking
Access Port	1/2" Diameter Top Back Wall
Interior Finish	Stainless Steel
Exterior Finish	Stainless Steel
Casters	Four locking
Shelves	Four Adjustable
Interior Lighting	LED light
Insulation	Polyurethane foam
Refrigerant:	R290 CFC free
Air Circulation	Forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	6.0 / 15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	420 lbs.
Warranty-(Domestic)	Two year parts and labor Compressors additional 3 years, no labor.



Specifications subject to change without notice