

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DHK4-49SD

49 Cubic Ft. (1387 liters) Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F
Two Swinging Solid Doors

Specifications:

| | |
|----------------------|----------------------------|
| * Capacity | 49 cubic ft. / 1387 liters |
| * Temperature Range | 2°C to 8°C / 36°F to 46°F |
| * Dimensions | |
| * Exterior w x d x h | 56.69" x 31.5" x 79.25" |
| * Interior w x d x h | 51.97" x 26.18" x 59.84" |

Features:

| | |
|----------------------|--|
| Temperature Control | Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay Battery back up |
| Access Port | 1/2" Diameter |
| Doors | Two hinged solid, locking |
| Interior Finish | White coated steel |
| Exterior Finish | White coated steel |
| Casters | Four, Locking |
| Shelves | Eight, Adjustable |
| Interior Lighting | LED light |
| Insulation | Polyurethane foam |
| Refrigerant: | R290A |
| Compressor | 1/5 hp |
| Air Circulation | Forced air |
| Defrost | Cycle, Auto |
| Voltage: Domestic | 115/60/1 |
| Running Amps/Breaker | 6.0/15 |
| Supply Plug | 5-15P NEMA UL listed |
| Shipping Weight | 450 lbs. |
| Warranty-(Domestic) | Two year parts and labor Compressors additional 3 years, no labor. |

