

SO-LOW LABORATORY REFRIGERATOR MODEL DHK4-72GD-SS

**72 Cubic Ft. (2095 liters) Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F
Three Glass Hinged Doors**

Specifications:

* Capacity	72 cubic ft. / 2095 liters
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	85.04" x 33.46" x 79.25"
* Interior w x d x h	80.31" x 26.18" x 59.84"

Features:

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay Battery back up
Access Port	1/2" Diameter Top Back Wall
Interior Finish	Stainless Steel
Exterior Finish	Stainless Steel
Doors	Three glass – hinged – key lock
Shelves	Twelve, adjustable
Casters	Four, Locking
Interior Lighting	LED light
Insulation	Polyurethane foam
Refrigerant	R290A CFC free
Compressor	2 x 1/3 hp
Air Circulation	Forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	12.0 /15
Supply Plug	NEMA 5-15P UL listed
Shipping Weight	700 lbs.
Warranty-(Domestic)	Two year parts and labor Compressors additional 3 years, no labor.

