

SO-LOW LABORATORY FREEZER TO -20°C MODEL DHN20-16MDF

**17.1 Cubic Ft. (484 liters) Upright – Temperature Range:
-20°C Solid Door**

Specifications:

* Capacity	17.1 cubic ft. / 484 liters
* Temperature Range	-10°C to -20°C
* Dimensions	
* Exterior w x d x h	32" x 30.5" x 64.5"
* Internal w x d x h	27.5" x 21.5" x 52"

Features:

- *MANUAL DEFROST
- *R134A REFRIGERATION SYSTEM
- *3 WHITE EVAPORATOR SHELVES
- *5 FULL- WIDTH FIXED DOOR BINS
- *1 FULL ACCESS SLIDING STORAGE BASKET
- *AUTOMATIC DOOR CLOSER
- * DEFROST DRAIN
- *INTERIOR LIGHT
- *ADJUSTABLE TEMPERATURE CONTROL
- *LOCK WITH POP OUT KEY
- *115 VOLT, 60 HERTZ, 1 PHASE
- *WARRANTIES: 12 MONTH PARTS

