

SO-LOW LAB & PHARMACY FREEZER MODEL DHN20-24MDS

"SETTING TEMPERATURE STANDARDS IN EVERY LAB"

24 Cubic Ft. Upright – Temperature Range: -10°C to -25°C

Specifications:

* Capacity	24 cubic ft.
* Temperature Range	-10°C to -25°C
* Dimensions	
* Exterior w x d x h	27.5" x 34.875" x 79.625"
	69 x 88 x 202 (cm)
Interior w x d x h	23.5" x 30" x 59"
	59 x 76 x 150 (cm)

Features:

Control	Mechanical
Doors	Internal chamber Thermometer
Interior finish	Single Solid, locking
Casters	White
Shelves	Four
Interior Lighting	Three
Insulation	Recessed
Refrigerant:	Polyurethane
Compressor	R404A CFC free
Condenser	3/4 hp
Evaporator	Fin and tube
Thermostat	Cold Wall
Defrost	Adjustable
Voltage: Domestic	Manual
Running Amps/Breaker	115/60/1
Supply Plug	12/15
	5-15P NEMA
	UL listed, CSA approved
Warranty-(Domestic)	1 ½ years parts and labor.
	Compressors additional
	3 ½ years, no labor.



Specifications subject to change without notice