

ENGINEERS AND MANUFACTURERS OF ULTRA-LOW FREEZERS

## ***SO-LOW FLAMMABLE MATERIAL STORAGE REFRIGERATOR MODEL DHN4-24FMS***

**24 Cubic Ft. Upright – Temperature Range: 2°C to 10°C**

**Solid Door**

### ***Specifications:***

|                      |   |
|----------------------|---|
| *Capacity            | 24 cubic ft.                                    |
| *Temperature Range   | 2°C to 10°C                                     |
| *Dimensions          |   |
| * Exterior w x d x h | 27.5" x 34.875" x 79.625"<br>69 x 88 x 202 (cm) |
| * Interior w x d x h | 23.5" x 30" x 59"<br>59 x 76 x 149 (cm)         |

***For Storage of flammables in approved containers  
for Class 1, Division 2 Storage of Flammable  
Liquids per NFPA Standards, 70, 99, 45, 497A, UL471***

### ***Features:***

|                      |  |
|----------------------|--|
| Controller           | Mechanical   |
| Doors                | Single Solid, locking  |
| Interior finish      | Smooth Aluminum  |
| Casters              | Four   |
| Shelves              | Three  |
| Insulation           | Polyurethane   |
| Refrigerant:         | R404A CFC free   |
| Condenser            | Fin and tube   |
| Evaporator           | Cold Wall  |
| Defrost              | Manual   |
| Voltage: Domestic    | 115/60/1   |
| Running Amps/Breaker | 6.4/15   |
| Supply Plug          | 5-15P NEMA<br>UL listed  |
| Warranty-(Domestic)  | 1 ½ years parts and labor.<br>Compressors additional<br>3 ½ years, no labor. |



Specifications subject to change without notice