

SO-LOW LAB & PHARMACY REFRIGERATOR WITH DRAWERS MODEL DHN4-24GD-D

“SETTING TEMPERATURE STANDARDS IN EVERY LAB”

24 Cubic Ft. Upright – Temperature Range: 2°C to 10°C / Glass Door

Specifications:

*Capacity	24 cubic ft.
*Temperature Range	2°C to 10°C
*Dimensions	
* Exterior w x d x h	27.5” x 35.5” x 79.625” 70 x 90 x 207 (cm)
* Interior w x d x h	23.5” x 30” x 59” 60 x 76 x 150 (cm)

Features:

Temperature Control	Digital display with Hi/Low alarm
Doors	Single glass, locking
Interior finish	White
Casters	Four
Drawers	4 stainless steel drawers
Interior Lighting	Fluorescent
Insulation	Polyurethane foam
Refrigerant:	R134A CFC free
Compressor	1/4 hp
Condenser	Fin and tube
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Thermostat	Adjustable
Defrost	Cycle, automatic
Voltage: Domestic	115/60/1
Running Amps/Breaker	8.6/15
Supply Plug	5-15P NEMA
Shipping Weight	365 lbs.
Warranty-(Domestic)	1 ½ years parts and labor. Compressors additional 3 ½ years, no labor.

