

SO-LOW LAB & PHARMACY REFRIGERATOR WITH DRAWERS MODEL DHN4-33GD-D

“SETTING TEMPERATURE STANDARDS IN EVERY LAB”

33 Cubic Ft. Upright – Temperature Range: 2°C to 10°C / Glass Door

Specifications:

* Capacity	33 cubic ft.
* Temperature Range	2°C to 10°C
* Dimensions	
* Exterior w x d x h	31.75" x 36.5" x 87.625" 80 x 92 x 222 (cm)
* Interior w x d x h	27.75" x 31" x 67" 70 x 79 x 170 (cm)

Features:

Temperature Control	Digital display with Hi/Low Alarm
Doors	Single glass, locking
Interior finish	White
Casters	Four
Drawers	Four Stainless Steel (18.5" x 22" x 3")
Interior Lighting	Fluorescent
Insulation	Polyurethane foam
Refrigerant:	R134A CFC free
Compressor	1/4 hp
Condenser	Fin and tube
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Thermostat	Adjustable
Defrost	Cycle, automatic
Voltage: Domestic	115/60/1
Running Amps/Breaker	8.6/15
Supply Plug	5-15P NEMA UL listed, CSA approved
Shipping Weight	439 lbs.
Warranty-(Domestic)	1 ½ years parts and labor. Compressors additional 3 ½ years, no labor.

