

SO-LOW UNDERCOUNTER REFRIGERATOR MODEL DHN4-5MDUCR

“SETTING TEMPERATURE STANDARDS IN EVERY LAB”

5.1 Cubic Ft. Upright – Temperature Range 2°C to 8°C

Specifications:

* Capacity	5.1 cubic ft
* Temperature Range	2°C to 8°C
* Dimensions	
* Exterior w x d x h	27" x 25" x 34" 69 x 64 x 86 (cm)
* Interior w x d x h	23" x 19.5" x 22.88" 58 x 50 x 58 (cm)

Features:

- *Manual Defrost
- *Air cooled condensing unit
- *1 Epoxy coated shelf
- *Door lock
- *Expansion valve system
- *White exterior, white poly interior
- *Gray epoxy coated shelf
- *Non-CFC refrigerant R134A
- *115 Volt, 60 Hertz, 1 Phase. UL & C-UL listed 3.5 Total amp draw. NEMA 5-15P plug provided.
- *Warranty, 18 months parts & labor, 5 years on compressor



Specifications subject to change without notice