

SO-LOW LAB & PHARMACY FREEZER MODEL DHN40-21SD

21 Cubic Ft. Upright – Temperature Range: -30°C to -40°C

Solid Door

Specifications:

* Capacity	21 cubic ft.
* Temperature Range	-30°C to -40°C
* Dimensions	
* Exterior w x d x h	27.5" x 35.75" x 81.5" 69 x 90 x 207 (cm)
* Interior w x d x h	23.5" x 30" x 59" 59 x 76 x 149 (cm)

Features:

Temperature Control	Digital display with Hi/Low alarm
Doors	Single solid, locking
Interior finish	White
Casters	Four
Shelves	Three
Insulation	Polyurethane foam
Refrigerant:	R404A
Compressor	1 hp
Condenser	Fin and tube
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Thermostat	Adjustable
Defrost	Cycle, automatic
Voltage: Domestic	115/60/1
Running Amps/Breaker	12.6/20
Supply Plug	5-20P
	UL/CUL listed
Shipping Weight	485 lbs.
Warranty-(Domestic)	1 ½ years parts and labor. Compressors additional 3 ½ years, no labor.

