

MODEL DHP4-66SGD

TOUCHSCREEN LAB REFRIGERATOR 2°C TO 8°C

Capacity [Cu. Ft.]	66 Cu. Ft.
Capacity [Liters]	1868 Liters
Temperature Range [°C]	2°C to 8°C
Temperature Range [°F]	36°F to 46°F
Dimensions	Width x Depth x Height
Exterior [Inches]	75" x 32" x 81"
Interior [Inches]	72" x 27" x 55"
Exterior [Centimeters]	190 x 81 x 206
Interior [Centimeters]	182 x 69 x 140
Compressor	Hermetic Compressor
Air Circulation	Forced Air
Defrost Method	Cycle
Refrigerant	R290
Exterior	Painted Steel
Interior	Painted Steel
Insulation	Polyurethane Foam
Lockable Door	3 Sliding Doors, Keyed Lock
Door Material	Glass
Interior Light	Adjustable LED
Porthole	1/2" Diameter
Casters	4, Locking
Storage	15 Adjustable Shelves
Shipping Weight	850 lbs.

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 6.2A
NEMA Configuration	NEMA 5-15
Certification(s)	UL Listed

Controller

Display Type	Microprocessor, Digital
Touchscreen	Yes, 4.3" Color Display
Alarms	High/Low Temp
Alarm Notification	Yes, Audible, Email, Text
Data Logging	Yes, Adjustable Intervals
Data Download	Yes, via USB, CSV Format
LAN Remote Access	Yes, via VNC
Min / Max Temperature	Yes, Display and Reset
Alarm Relay	Yes, Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 7 Years Parts, Labor & Compressor Part Warranty.