

## *SO-LOW LAB & PHARMACY FREEZER*

### *MODEL DHS25-25GD*

*“SETTING TEMPERATURE STANDARDS IN EVERY LAB”*

**25 Cubic Ft. (679 liters) Glass Door Upright**

**Temperature Range: -18°C to -25°C / 0°F to -13°F**

#### **Specifications:**

* Capacity	25 cubic ft. / 679 liters
* Temperature Range	-18°C to -25°C / 0°F to -13° F
* Dimensions	
* Exterior w x d x h	28.35” x 31.5” x 79.25”
* Interior w x d x h	23.62” x 26.18” x 59.84”

#### **Features:**

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay Battery back up
Access Port	1/2" Diameter
Door	Single glass – lockable
Casters	heavy duty, swivel
Shelves	3 epoxy coated adjustable
Interior Lights	Led light
Insulation	Polyurethane CFC Free
Refrigerant:	R-290
Air Circulation	Forced air
Evaporator	Fin and tube, forced air
Thermostat	Adjustable
Defrost	Automatic
Voltage:	115/60/1
Running Amps/Breaker	10.3/15
Supply Plug	5-15 NEMA
Shipping Weight	255 lbs.
Warranty-(Domestic)	Two year parts and labor Compressors additional 3 years, no labor.



Specifications subject to change without notice.