

## ***SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DHS4-49GD-SS***

**49 Cubic Ft. (1387 liters) Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F**

**Two Swinging Glass Doors**

### ***Specifications:***

* Capacity	49 cubic ft. / 1387 liters
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	56.69" x 31.5" x 79.25"
* Interior w x d x h	51.97" x 26.18" x 59.84"

### ***Features:***

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay Battery back up
Access Port	1/2" Diameter Top Back Wall
Doors	Two hinged glass, locking
Interior Finish	Stainless Steel
Exterior Finish	Stainless Steel
Casters	Six, locking
Shelves	Eight, Adjustable
Interior Lighting	LED light
Insulation	Polyurethane foam
Refrigerant:	R290A
Compressor	1/5 hp
Air Circulation	Forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	6.0/15
Supply Plug	5-15P NEMA
	UL listed
Shipping Weight	450 lbs.
Warranty-(Domestic)	Two year parts and labor Compressors additional 3 years, no labor.



Specifications subject to change without notice