

## *SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DHS4-49SGD-SS*

**49 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F**

**Two Siding Glass Doors**

### *Specifications:*

* Capacity	49 cubic ft. / 1286 liters
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	56.69" x 31.5" x 79.25"
* Interior w x d x h	51.97" x 26.18" x 59.84"

### *Features:*

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay Battery back up
Access Port	1/2" Diameter Top Back Wall
Interior Finish	Stainless Steel
Exterior Finish	Stainless Steel
Doors	Two sliding glass
Casters	Six, locking
Shelves	Eight, Adjustable
Interior Lighting	LED light
Insulation	Polyurethane foam
Refrigerant:	R290A
Compressor	1/5 hp
Air Circulation	Forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	6.0/15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	450 lbs.
Warranty-(Domestic)	Two year parts and labor Compressors additional 3 years, no labor.

