

SO-LOW LAB & PHARMACY REFRIGERATOR MODEL DHS4-49SGD

49 Cubic Ft. Upright – Temperature Range: 2°C to 8°C / 36°F to 46°F

Two Siding Glass Doors

Specifications:

* Capacity	49 cubic ft. / 1286 liters
* Temperature Range	2°C to 8°C / 36°F to 46°F
* Dimensions	
* Exterior w x d x h	56.69" x 31.5" x 79.25"
* Interior w x d x h	51.97" x 26.18" x 59.84"

Features:

Temperature Control	Digital Display With Hi/Low Alarm Alarm Relay Dry Contacts Min/Max Memory Power Failure Relay Battery back up
Access Port	1/2" Diameter Top Back Wall
Interior Finish	White coated steel
Exterior Finish	White coated steel
Doors	Two sliding glass
Casters	Six, locking
Shelves	Eight, Adjustable
Interior Lighting	LED light
Insulation	Polyurethane foam
Refrigerant:	R290A
Compressor	1/5 hp
Air Circulation	Forced air
Defrost	Cycle, Auto
Voltage: Domestic	115/60/1
Running Amps/Breaker	6.0/15
Supply Plug	5-15P NEMA UL listed
Shipping Weight	450 lbs.
Warranty-(Domestic)	Two year parts and labor Compressors additional 3 years, no labor.

