

SO-LOW LABORATORY FREEZER TO -20°C PREMIER MODEL DHW20-20MDP

19.7 Cubic Ft.(558 liters) Upright

Temperature Range: -20°C

Solid Door

Specifications:

* Capacity		19.7 cubic ft. / 558 liters
*Temperature Range		-20°C
* Dimensions		
* Exterior	w x d x h	32" x 28.75" x 82"
* Internal	w x d x h	27" x 20" x 57"

Features:

CONTROL

Digital temperature control with two temperatures displays. One display shows the set temperature, the other shows the actual temperature in the freezer.

ALARM

The battery operated alarm system will emit an audible and visual signal when there is a mechanical or electrical failure. The alarm has an over and under temperature setting, alarm silencing switch, and battery test switch.

REMOTE ALARM RELAY

A relay for a remote alarm hook-up is also provided.

Manual Defrost

R134A Refrigeration System

6 In-Door Shelves

4 Fixed Shelves

Bottom Gliding Basket

Interior Light

Key Lock

115 Volt, 60 Hertz, 1 Phase

One Year Warranty

Optional 7-day Temperature Recorder

With 6" Charts

