

MODEL DHW20-30MDP

LAB AND PHARMACY FREEZER 0°C TO -20°C

Capacity [Cu. Ft.]	30 Cu. Ft.
Capacity [Liters]	849 Liters
Temperature Range [°C]	0°C to -20°C
Temperature Range [°F]	32°F to -4°F

Dimensions	Width x Depth x Height
Exterior [Inches]	35.5" x 35" x 80"
Interior [Inches]	32" x 22.5" x 65"
Exterior [Centimeters]	90 x 88 x 203
Interior [Centimeters]	81 x 57 x 165

Compressor	Hermetic Compressor
Air Circulation	Gravity Flow
Defrost Method	Manual
Refrigerant	R404A
Exterior	Painted Steel
Interior	Painted Steel
Insulation	CFC-Free
Lockable Door	Keyed
Porthole	Yes, 3/8" Diameter
Door Handle	Stainless Steel
Leveling Legs	4
Storage	4 Fixed Shelves

Electrical

Voltage:	115 Volts 60 Hertz 1 Phase
Amperage Line Running	15 Amp Dedicated 6A
NEMA Configuration	NEMA 5-15

Controller

Display Type	Microprocessor, Digital
Alarms	High/Low Temperature
Min / Max Temperature	Yes, Display and Reset
Alarm Relay	Yes, Dry Contacts

IMAGE BELOW DOES NOT SHOW OPTIONAL ACCESSORIES



Warranty 1 Year Parts & Labor